INSTRUCTIONS FOR USE

ISOLITE[®]

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TR

Device Description

Isolite® Illuminated Dental Isolation System

REF EIL0201 (USA)

REF EIL0204 (EU)

REF EIL0205 (AUS)

This system is intended for use with the Mouthpieces developed by Zyris, Inc. These Mouthpieces are available in a wide range of sizes for effective isolation.

Each Isolite Illuminated Dental Isolation System is comprised of four main components: 1) Titanium Control Base, 2) Vacuum Light Pipe, 3) Power Hose, 4) LED Smart Stick.

Titanium Control Base

REF PIL0201

The Titanium Control Base includes two Vacuum Levers (P/N: PIL0205) to adjust suction in the upper and lower areas of the Mouthpiece.

LED Smart Stick

REF PIL0202

The LED Smart Stick controls the LED light source. The one-button light control adjusts: on/off, 5 levels of brightness, and a cure-safe mode.

Vacuum Light Pipe

REF PIL0204

The Vacuum Light Pipe transmits ultra-bright light while channeling suction. The Mouthpiece attaches onto the Vacuum Light Pipe.

Power Hose

REF PIL0203 (USA) (1.82 m / 72 in)

REF PIL0208 (USA) (0.91 m / 36 in)

REF PIL0209 (USA) (1.22 m / 48 in)

REF PIL0210 (USA) (1.52 m / 60 in)

REF PIL0211 (EU/AUS) (1.32 m / 52 in)

The Power Hose delivers both power and vacuum suction from the vacuum canister to the Control Head. The Power Hose connects to the vacuum canister. At the end of the Power Hose, an indicator light will alert the user if an incorrect power adapter has been plugged into the Power Hose.

Power Adapter

REF PIL0207

The Power Adapter delivers electrical power to the Power Hose and the LED Smart Stick. It connects to the Power Hose at the Power Jack and plugs into a standard power outlet (110v or 220v).

Plug Adapter kit

REF AIL0219 (EU)

REF AIL0220 (NZ/AUS)

International power plug adapters.

Suction Line Plug

REF PIL0206

O-ring Maintenance Kit

REF AIL0206 (USA/AUS)

REF AIL0223 (EU)

3/8" Y-Adapter Kit

REF AIL0207 (USA/AUS)

9.5 mm Y-Adapter Kit

REF AIL0222 (EU)

Double End Tube Brush

REF AIL0204

Mouthpiece Lubricant

REF AIL0202

Barrier Sleeves

REF AIL0203 (non-sterile)

Sold Separately

Non-sterile Control Head Barrier Sleeves are available for additional protection from cross-contamination. Use of Barrier Sleeves does not eliminate the need to follow disinfecting and sterilization recommendations for the Control Head and Vacuum Light Pipe.

Mouthpieces

REF CILXX0X (non-sterile)

Sold Separately

Mouthpieces are single use only and available in various sizes. The Mouthpieces are constructed of a soft, flexible polymer and are not made with natural rubber latex. They do not contain phthalates or BPA's. For more information on sizes, go to Zyris.com.

Intended Use

Isolite is a dental isolation system for clinical use by dental professionals. The system is designed to isolate two quadrants at once, simultaneously providing intra-oral illumination, tongue and cheek retraction, and continuous suction. The retraction and continuous suction are the important functions that control the oral environment, reducing contamination and improving bond strength while protecting the airway and soft tissue.

Contraindications

No known contraindications. Do not use without reviewing how to select the correct Mouthpiece size for each patient. Do not use if unable to position the Mouthpiece correctly.

Side Effects

Use of the Mouthpiece may trigger a pharyngeal reflex (laryngeal spasm) in hypersensitive patients, or if it is not placed correctly within the oral cavity.

How Supplied

Isolite

The Isolite Illuminated Dental Isolation System contains: Titanium Control Base, LED Smart Stick, 6 Vacuum Light Pipes, Power Hose, Power Supply, Plug Adapter Kit (for specific countries), Instructions for Use.

Accessories: Suction Line Plug, Mouthpiece Lubricant, Y-adapter Kit, Double End Tube Brush, and O-ring Maintenance Kit.

Mouthpieces

Mouthpieces (sold separately) are supplied in convenient operatory boxes. Each Mouthpiece is individually packaged.

Installation

Isolite is designed to install easily in most dental operatories (Figure 1). Installation and usage instructions, including a video demonstration, are also available at zyris.com.

Step 1 Connect the Power Hose to the vacuum canister by fitting the end of the vacuum hose onto an auxiliary HVE port on the vacuum canister. If an auxiliary HVE port is not available, use the provided Y-adapter, so that the standard HVE hose may remain connected.


Figure 1. Installation

NOTE: For countries outside the USA, installation may need to be completed by a qualified service technician. Please consult your in-country representative for additional information.

Step 2 Plug the Power Adapter into an electrical outlet. Plug the DC connector into the Power Jack at the end of the Power Hose. Connect the Control Head to the head of the Power Hose. When properly installed, a green LED indicator light will appear on the LED Smart Stick. Once connected, test the dual Vacuum Levers to confirm suction functionality.

NOTE: For countries outside the USA, attach the appropriate plug adapter to the power adapter. See the illustration included with the Plug Adapter Kit. Direct wire installations require completion by a qualified service technician. Please consult your in-country representative for additional information.

Step 3 Place the Isolite in an open slot on the toolbar when not in use.

How to Operate

Using the Dual Vacuum Levers

Dual Vacuum Levers adjust suction intensity in both the upper and lower quadrants. To adjust, rotate the left and right Vacuum Levers to the desired suction intensity for each quadrant. To shut off suction, rotate the levers to the most downward position.

Using The Led Smart Stick

Adjust lighting functions by using the one-button light control.

Operation	Action
Turn light on / off	Tap button once
Set brightness	Press and hold button to cycle through brightness levels; release to set
Engage / disengage cure-safe mode	With light on, double-tap button

Table 1. LED Operation

The LED Smart Stick features an LED indicator which provides the status of the system.

LED Indicator	Status
Steady green, pulsing every 10 sec	System is ready; power is connected and the light is off
Steady green	Light is on
Blinking blue	With button pressed and held, marks each of the 5 brightness levels
Steady amber	Cure-safe mode
Steady blue	Self-regulating temperature control has automatically lowered the brightness to reduce heat level

Table 2. LED Status

Selecting the Mouthpiece

To select the correct Mouthpiece size, measure the patient's interincisal opening. If between sizes, you may select the smaller size.

Mouthpiece Size	Approx. Interincisal Opening
Pediatric (PED)	< 30 mm
Extra Small (XSM)	< 30 mm
Small (SML)	30 mm
Medium (MED)	30-45 mm
Medium Deep Vestibule (MDV)	30-45 mm
Large (LRG)	> 45 mm

Table 3. Mouthpiece Sizes

The Finger Method

The finger method may also be used to select the appropriate Mouthpiece size (Figure 2). For additional assistance, watch a video demonstration online at Zyris.com.


Figure 2. The Finger Method



Figure 3. Mouthpiece

Attaching the Mouthpiece

Slide the Mouthpiece onto the Vacuum Light Pipe until the two locking barbs engage. Water or mouthpiece lubricant may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.

Placing the Mouthpiece

Prior to placing the Mouthpiece within the patient's oral cavity, moisten the Mouthpiece and the patient's lips.

Step 1

Step 2

Step 3

Step 1

Hold the sides of the Control Head with your index finger and thumb. With your other hand, fold the Cheek Shield onto the Tongue Retractor. Gently slide the folded mouthpiece into the buccal vestibule on the side to receive treatment. Angle the mouthpiece to allow the lower edge of the Tongue Retractor to move along the buccal edge of the teeth.

Step 2

Move the bite block onto the occlusal surface of the teeth, just distal to the mandibular cuspid. Instruct the patient to "rest gently" on the bite block to secure.

Step 3

Place the Cheek Shield into the buccal vestibule. Move the Isthmus onto the retromolar pad behind the maxillary tuberosity. Adjust the Tongue Retractor in the lingual vestibule as needed. Move the bite block distally to provide more vertical working room.

Figure 4, 5, & 6. Mouthpiece Placement

Disassembly of the Control Head

Step 1

Press the release button (Figure 7) on the underside of the Titanium Control Base to release and disconnect the Vacuum Light Pipe.

Step 2

Remove the LED Smart Stick from the Titanium Control Base.

Step 3

Holding the head of the Power Hose, pull to remove the Titanium Control Base. Do not twist.


Figure 7. Control Head

Disinfection, Cleaning, and Sterilization

Surface Disinfection of the Titanium Control Base

Zyris, Inc. recommends performing a surface disinfection of the Titanium Control Base between each patient. To disinfect the Titanium Control Base, first remove the Vacuum Light Pipe and the LED Smart Stick.

If the Titanium Control Base is visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Otherwise, use a hospital grade disinfectant and follow the manufacturer's instructions for contact time.

Surface Disinfection of the LED Smart Stick

Zyris, Inc. recommends performing a surface disinfection of the LED Smart Stick between each patient.

Spray a ready-to-use, fast-killing disinfectant onto entire surface of the LED Smart Stick and allow to remain wet per manufacturer's instructions. Using a nylon bristle brush, clean around the light control button area. Do not press the button while cleaning. Gently work the bristles in a circular motion, following the groove between the housing and the button for approximately 60 seconds to remove all visual debris. Repeat if any contaminants are still visible.

Wipe the LED Smart Stick thoroughly with a ready-to-use, fast-killing disinfectant for the specified contact time per manufacturer's instructions. Dry off the LED Smart Stick and gold connector pins before use.

Surface Disinfection of the Power Hose

Zyris, Inc. recommends performing a surface disinfection of the Power Hose between each patient.

Use a hospital grade disinfectant per manufacturer's instructions. If components are visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Barrier Sleeves may also be used and changed between patients.

Cleaning and Sterilization of the Vacuum Light Pipe

Zyris, Inc. recommends sterilization of the Vacuum Light Pipe between each patient.

Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not use solution to hold parts. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 4).

Cleaning and Sterilization of the Titanium Control Base

Zyris, Inc. recommends following normal infection control protocol for dental devices. Users may sterilize the Titanium Control Base at their discretion.

Remove the Vacuum Levers from the Titanium Control Base. Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not hold parts in cleaning solution or use ultrasonic cleaners. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 4).

After sterilization, lubricate the Vacuum Lever O-rings with the Dow 111 O-ring lubricant provided in the O-ring Maintenance Kit.

Validated Sterilization Cycles

Cycle	Temperature	Min. Exposure Time	Min. Drying Time
Standard	121 °C / 250 °F	30 min	30 min
Standard	132 °C / 270 °F	15 min	15 min
Prevacuum	132 °C / 270 °F	4 min	20 min
Prevacuum	134 °C / 273 °F	3 min	20 min

Table 4. Validated Sterilization Cycles

Maintenance

Titanium Control Base

The Vacuum Lever O-rings require weekly lubrication to operate at peak efficiency. Remove the Titanium Control Base from the Power Hose. To remove the Vacuum Levers, push the levers to the “off” position and keep pushing gently until each lever pops out, or gently pry them out. Coat the O-ring by applying a small amount of the Dow 111 O-ring Lubricant provided in the O-ring Maintenance Kit. Apply an additional amount to the Titanium Control Base retainer rings as well. Reattach the Vacuum Levers to the Titanium Control Base (Figure 8 and 9).

Figure 8. Vacuum Levers

Figure 9. Titanium Control Base

Power Hose

Clean the inside of the Power Hose regularly according to frequency of use to prevent the accumulation of biofilm. Start by vacuuming water through the Power Hose. Then vacuum a dental hose cleaner through the hose per manufacturer's instructions. The vacuum hose cleaner should be non-foaming, deodorizing, anti-microbial, and able to dissolve organic debris. Most enzymatic cleaners with a pH between 6 and 8 are acceptable. Do not use disinfectants as a vacuum hose cleaner. If further disinfection is necessary, consult your vacuum manufacturer for recommendations.

To lubricate the Power Hose O-ring, remove the Titanium Control Base from the Power Hose. The black O-ring is located approximately 13 mm (½ inch) down the inside of the blue fitting at the head of the Power Hose. (Figure 10) Swab the O-ring with the Dow 111 O-ring lubricant provided in the O-ring Maintenance Kit.

Figure 10. Power Hose

Warnings

⚠ WARNING: Shock hazard if incorrect power supply is used. Only use the power adapter supplied with your Isolite. Use of other power sources may damage electronics and voids your product warranty. The electrical system requirement is 5.0 VDC, 1.0 amp.

⚠ WARNING: Do not place a metal HVE valve onto the Power Hose. This may damage electrical components and may present a shock hazard.

⊗ WARNING: Mouthpieces are single use only. Reuse can result in the transmission of infectious pathogens. Mouthpieces must be discarded after use. Mouthpieces cannot be safely resterilized after use.

⊗ WARNING: Do not use if the individual Mouthpiece package has been previously damaged or opened. Mouthpieces are individually packaged in polyethylene bags.

⚠ WARNING: Be sure to use established protocols and patient monitoring when using the Mouthpiece in conjunction with sedation dentistry. Patient management by a certified anesthesiologist is strongly recommended.

Precautions

⚠ CAUTION: When performing dental procedures requiring high levels of blood evacuation, it is recommended that the Mouthpiece be intermittently flushed with water to prevent coagulation in the Mouthpiece to avoid loss of suction and light transmission.

⚠ CAUTION: Do not block access to the Power Adapter. The Power Adapter should be easily accessible so that it can be unplugged from the wall outlet if necessary.

⚠ CAUTION: To avoid electrical shock, the use of an extension cord with the Power Adapter is not recommended.

⚠ CAUTION: Do not crimp or fold the Power Hose to shut off suction. Internal damage to the Power Hose wiring may occur.

⚠ CAUTION: The LED Smart Stick contains sensitive electronics and requires special care and handling. Do not submerge, autoclave, Chemclave, or place in ultrasonic cleaner.

⚠ CAUTION: The LED Smart Stick should only be removed and installed while the Titanium Control Base is attached to the power hose. Failure to do so may result in damage to the gold connector pins.

⚠ CAUTION: Do not autoclave the Titanium Control Base at temperatures above 135 °C; higher temperatures could damage components and void your product warranty.

⚠ CAUTION: Do not press the LED Light Control Button while cleaning.

⚠ CAUTION: Do not hold Titanium Control Base or LED Smart Stick in cleaning solution or use ultrasonic cleaners.

⚠ CAUTION: Failure to maintain properly lubricated O-rings can result in damage to the Isolite that is not covered by product warranty.

⚠ CAUTION: With regular use, cleaning, and maintenance, the device is designed to last for years. Check the device for wear prior to use to make sure that it is adequate for use.

Guide to Symbols

Non-sterile

Catalog number

Caution

Do not re-use

Positive polarity

Consult Instructions for Use

Do not use if package is damaged


Type B applied part electrical hazard

Date of manufacture

Manufacturer

Storage Conditions

Limited 1-Year Warranty

Zyris Inc. warrants this product to be free from defects in material and workmanship for a period of one year from the date of original purchase. If the product exhibits such a defect, Zyris, Inc. will, at its option, repair or replace it without costs for parts and labor. The product must be returned by the customer to Zyris, Inc. in accordance with current Return Material Authorization procedures. Proof of date of original purchase may be required.

This warranty does not cover finishes or normal wear, nor does it cover damage resulting from accident, misuse, dirt, tampering, unreasonable use, failure to provide reasonable and necessary maintenance, service performed or attempted by unauthorized individuals or service agencies, or devices that have been otherwise modified. All implied warranties, including any implied warranty or merchantability or fitness for any particular purpose, are limited in duration to one year from the date of receipt of the product. In no event will Zyris, Inc. be responsible for consequential damages resulting from the use of this product.

Instructions for Use may also be found at:
Zyris.com

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All rights reserved.